## REFERENCE VALUE

ACL TOP Family/ACL TOP Family 50 Series<sup>†</sup>

ng/mL

C

LOT

C

## Hemosil®

## D-Dimer HS 500 - 0020500100

## Intended use

Hemosil D-Dimer HS 500 is an automated latex enhanced immunoassay for the quantitative determination of D-Dimer in human citrated plasma on the ACL TOP Family and ACL TOP Family 50 Series<sup>†</sup> Systems for use, in conjunction with a clinical pretest probability (PTP) assessment model to exclude venous thromboembolism (VTE) in outpatients suspected of deep venous thrombosis (DVT) and pulmonary embolism (PE).

## Summary and principle

D-Dimer is contained in the soluble derivative formed upon plasmin degradation of Factor XIIa cross-linked fibrin (FDP). Plasmin, a serine protease, when free from inhibitors digests the insoluble cross-linked fibrin yielding a variety of soluble derivatives. One of these derivatives is the Fibrinopeptide fragment (FDP), which contains the D-Dimer domain. These soluble fragments contain a neopeptide (D-Dimer domain) which is not present on the original fibrinogen molecule, its degradation products or on soluble fibrin.<sup>1,2</sup> Determination of D-Dimer is becoming a widespread tool for diagnosing thrombosis and monitoring anticoagulation therapy.

Elevated levels of D-Dimer are seen in clinical conditions such as deep vein thrombosis (DVT), pulmonary embolism (PE) and during normal intramammary coagulation (DIC).

D-Dimer also rises with age<sup>3</sup> and normal pregnancy, but very high levels are associated with complications.<sup>4</sup>

Further, age-adjusted cut-off values for DVT and PE suspicion have been shown to increase specificity of D-Dimer and reduce the number of unnecessary imaging studies in patient populations greater than 50 years.<sup>5,6,7</sup>

A negative D-Dimer result when combined with a clinical assessment of pretest probability has been shown to have a high negative predictive value (NPV) for DVT and PE.<sup>8,9</sup>

The Latex Reagent is a suspension of polystyrene latex particles of uniform size coated with the F(ab)<sub>2</sub> fragment of a monoclonal antibody highly specific for the D-Dimer domain included in fibrin soluble derivatives. The use of the F(ab)<sub>2</sub> fragment allows a more specific D-Dimer detection avoiding the interference of some endogenous factors like the Rheumatoid Factor. When a plasma containing D-Dimer is mixed with the Latex Reagent and the Reaction Buffer included in the kit, the D-Dimer will bind to the latex particles, causing agglutination. The degree of agglutination is directly proportional to the concentration of D-Dimer in the sample and is determined by measuring the decrease of transmitted light caused by the aggregates (luminescence immunoassay).<sup>10</sup>

**Composition**

The D-Dimer HS 500 kit consists of:

**R Latex Reagent** (Cat. No. 0020500110): 3 vials x 4 mL of a suspension of polystyrene latex particles coated with the F(ab)<sub>2</sub> fragment of a mouse monoclonal antibody (MA-803) directed against D-Dimer containing bovine serum albumin, buffer, stabilizers and preservative.

**B Reaction Buffer** (Cat. No. 0020500120): 3 vials x 6 mL of Hepes buffer containing bovine serum albumin, buffer, stabilizers and preservative.

**C D-Dimer Calibrator** (Cat. No. 0020500130): 2 vials x 1 mL of a lyophilized solution of D-Dimer partially purified from human fibrin digested with human plasmin containing bovine serum albumin, buffer, stabilizers and preservative.

**PRECAUTIONS AND WARNINGS**

The following reagents are supplied by FDA approved test methods and found nonreactive for Hepatitis B Surface Antigen (HBsAg), Anti-HCV and HIV 1/2 antibodies. Handle as if potentially infectious.<sup>11</sup>

The Reaction Buffer contains less than 0.1% sodium azide that may form explosive azides in metal plumbing. Use proper disposal procedures.

**Hazard class:** none

**Hazard statement:** none

**Precautory statements:** none

**Latex Reagent:** 2.3% of the mixture consists of component of unknown acute toxicity (inhalation) for the human health.

**Reaction Buffer:** Up to 5.7% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

**D-Dimer Calibrator:** Up to 5.6% of the mixture consists of component of unknown acute toxicity (dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

This product is for *in vitro* Diagnostic Use.

**Precautionary statements:** none

**Latex Reagent:** Gently swirl several times to mix before use. Do not shake. Avoid foam formation.

**Reaction Buffer:** Gently mix several times before use. Do not shake. Avoid foam formation.

**D-Dimer Calibrator:** Dissolve the contents of each vial with 1 mL of CLSI (formerly NCCLS) CLRW Type water or equivalent.<sup>12</sup> Replace the stopper and gently mix. Make up of the complete reconstitution of the product. Keep the reagent at 15–25°C for 30 minutes and gently invert to mix before use. Do not shake. Avoid foam formation.

**Note:** Bubbles on top of the liquid may interfere with the instruments liquid sensors.

**Reagents storage and stability**

Unopened reagents and controls are stable until the expiration date shown on the vial when stored at 2–8°C.

**Latex Reagent:** Opened reagent is stable 1 month at 2–8°C in the original vial or 7 days at 15°C on the ACL TOP Family and ACL TOP Family 50 Series<sup>†</sup>.

**Reaction Buffer:** Opened reagent is stable 1 month at 2–8°C in the original vial or 7 days at 15°C on the ACL TOP Family and ACL TOP Family 50 Series<sup>†</sup>.

**PLEASE NOTE:** The D-Dimer HS 500 Latex Reagent and Reaction Buffer have the same stability claims and should be used as a pair in the same system. If one is replaced the other component (i.e. Latex Reagent or Reaction Buffer) should also be replaced at the same time regardless of the residual content in the vial.

**D-Dimer Calibrator:** Stability after reconstitution: 1 month at 2–8°C, 3 days at 15–25°C, or 2 months at -20°C in the original vial. Freezer Calibrator may be thawed at 37°C and gently mixed before use. Do not refreeze.

For optimal stability remove reagents from the system and store them at 2–8°C in the original vial.

⚠ Low relative humidity is associated with increased evaporation of uncapped reagents, which may decrease onboard stability. For optimum on-board stability, laboratory temperature and humidity should be controlled.

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