Free Protein S - 0020002700

Intended use

Automated latex ligand immunoassay for the quantitative determination of free Protein S in human citrated plasma on IL Coagulation Systems.

Summary and principle

Protein S is a vitamin K-dependent cofactor for the anticoagulant and the profibrinolytic effects of activated Protein C. Two forms of Protein S are present in plasma: free Protein S (40%), and Protein S linked to the complement C4b-binding protein (C4BP) (60%). Only free Protein S has functional cofactor activity.¹

Protein S deficiency may be hereditary¹ or acquired. Acquired deficiency may be observed during pregnancy, oral anticoagulant therapy, oral contraceptive use, in liver disease, in newborn infants as well as in other clinical conditions.^{4,5,6} Deficiency of Protein S has been associated with a high risk of developing venous thromboembolism especially in young people.^{7,8}

The presence of free Protein S is determined by measuring the increase of turbidity produced by the agglutination of two latex reagents. Purified C4BP adsorbed onto the first latex reagent reacts with a high affinity for free Protein S of patient plasma in the presence of Ca²⁺ ions.^{9,10,11} The free Protein S adsorbed on the C4BP latex triggers the agglutination reaction with the second latex reagent which is sensitized with a monoclonal antibody directed against human Protein S. The degree of agglutination will be directly proportional to the free Protein S concentration in the test sample.

Composition

The **Free Protein S** kit consists of:

B	C4BP Buffer (Cat. No. 0020002720): 3 vials x 4 mL of Borax buffer containing bovine serum albumin, stabilizers and preservative.		
R1	C4BP Latex (Cat. No. 0020002710): 3 vials x 4 mL of a lyophilized suspension of polystyrene latex particles coated with purified human C4BP containing bovine serum albumin, stabilizers and preservative.		
R2	Anti PS Mab Latex (Cat. No. 0020002730): 3 vials x 2 mL of a suspension of polystyrene latex particles coated with a monoclonal antibody directed against human Protein S containing bovine serum albumin, stabilizers and preservative.		

PRECAUTIONS AND WARNINGS:

The material in this product was tested by FDA approved test methods and found nonreactive for Hepatitis B Surface Antigen (HBsAg), Anti-HCV and HIV 1/2 antibodies. Handle as if potentially infectious.¹²

Danger

C4BP Buffer, C4BP Latex, ANTI PS Mab Latex

Hazard class: Reproductive Toxicity, cat 1B

Hazard statements: May damage fertility or the unborn child.

Precautionary statements: Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Wear protective gloves/protective clothing/eye protection/face protection. If exposed or concerned: Get medical advice/attention. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations.

Supplemental hazard information: EUH210: Safety data sheet available on request. This product is *For in vitro* Diagnostic Use.

Preparation

C4BP Buffer: The reagent is ready for use.

C4BP Latex: Dissolve the contents of each vial of C4BP latex by pouring the entire contents of one vial of C4BP buffer into the latex reagent vial. Replace the stopper and swirl gently for a minimum of 20 seconds to completely dissolve the lyophilized latex. Make sure of the complete reconstitution of the product. It must appear as a homogenous and slightly milky suspension. Keep the reagent at 15-25°C for 30 minutes and do not mix before use. Do not shake.

Anti PS Mab Latex: Invert to mix before use. Do not shake.

Note: Avoid foam formation when homogenizing reconstituted reagents. Bubbles on top of the liquids may interfere with the instruments liquid sensors.

Reagent storage and stability

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.

C4BP Latex - Stability after reconstitution: 1 month at 2-8°C in the original vial, 1 week at 15°C on the ACL Elite/Elite Pro and ACL TOP® Family and ACL TOP Family 50 Series¹. Do not freeze.

Anti PS Mab Latex - Opened reagent is stable 1 month at 2-8°C in the original vial, 1 week at 15-25°C on the ACL Elite/Elite Pro and 1 week at 15°C on the ACL TOP Family and ACL TOP Family 50 Series¹. Do not freeze.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

Low relative humidity is associated with increased evaporation of uncapped reagents, which may decrease onboard stability. For optimum on-board stability, laboratory temperature and humidity should be controlled.

Instrument/test procedures

Refer to the appropriate IL instrument's Operator's Manual and/or Application Manual for the complete assay procedure instructions.

Specimen collection and preparation

Nine parts of freshly drawn venous blood are collected into one part 3.2% trisodium citrate. Refer to CLSI Document H21-A5 for further instructions on specimen collection, handling and storage.¹³

Frozen plasma samples should be rapidly thawed at 37°C while gently mixing before testing. After thawing the assay must be performed within 2 hours.

Additional reagents and control plasmas

The following are not supplied with the kit and must be purchased separately.

	Americas and Pacific Rim	Europe
	Cat. No.	Cat. No.
Calibration plasma	0020003700	0020003700
Normal Control	0020003120/0020003110	0020003110
Special Test Control Level 1	0020011000	0020011000
Special Test Control Level 2	0020012000	0020012000
Low Abnormal Control	0020003220/0020003210	0020003210
Factor diluent	0009757600	0009757600

Quality Control

Normal and abnormal controls are recommended for a complete quality control program.¹⁴ Normal Control, Low Abnormal Control, and Special Test Control Levels 1 & 2 are designed for this program. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance to good laboratory practice. Refer to the instruments' Operator's Manual for additional information. Refer to Westgard *et al* for identification and resolution of out-of-control situations.¹⁵

Results

Free Protein S results are reported in % normality. Refer to the instrument's Operator's Manual for additional information. The assay results should be used with other information, including the clinical context, in forming a diagnosis.

Limitations/interfering substances

Free Protein S results on the ACL Systems are not affected by heparin (UF heparin or LMW heparin) up to 1.5 IU/mL, bilirubin up to 18 mg/dl, hemoglobin up to 200 mg/dl, lipids up to 1280 mg/dl, platelets up to 10¹¹/L and rheumatoid factor up to 350 IU/mL. Free Protein S results on the ACL TOP Family and ACL TOP Family 50 Series¹ are not affected by heparin (UF heparin or LMW heparin) up to 1.5 IU/mL, bilirubin up to 25 mg/dl, hemoglobin up to 200 mg/dl, triglycerides up to 1500 mg/dl, platelets up to 10¹¹/L and rheumatoid factor up to 900 IU/mL.

Free Protein S assay is not affected by Factor V Leiden mutation (APC-R).

ENGLISH - Insert revision 04/2019

Expected values

A normal range study was performed using the Free Protein S kit on IL Coagulation Systems.

Sex	N	ACL Elite/Elite Pro	N	ACL TOP Family	ACL TOP Family 50 Series ¹
Male (% free PS)	130	72.2 - 123.3	128	74.1 - 146.1	
*Female (% free PS)	102	57.6 - 112.5	128	54.7 - 123.7	
* Note: Age and hormonal status may affect the normal range for females. ¹⁶ Ranges were calculated as the mean of the two forms: Proteína S Libre (40%) y Proteína S ligada a la proteína transportadora de la fracción C4b del complemento (C4BP) (60%). Sólo la Proteína S libre tiene actividad funcional como cofactor. ⁷					
La deficiencia en Proteína S puede ser hereditaria ¹ o adquirida. La deficiencia adquirida se puede observar durante el embarazo, terapia anticoagulante oral, uso de anticonceptivos orales, en enfermedades hepáticas, en neonatos, así como en otras situaciones clínicas. ^{4,5,6} La deficiencia en Proteína S se ha asociado a un alto riesgo de desarrollar tromboembolismo venoso, especialmente en gente joven. ^{7,8}					

Performance characteristics

Precision:

Within run and total (run to run and day to day) precision was assessed over multiple runs.

ACL Elite/Elite Pro	Mean (% free PS)	CV% (Within run)	CV% (Total)
Normal Control	97.1	1.8	3.0
Special Test Control Level 1	56.6	2.9	4.9
Special Test Control Level 2	21.8	3.2	5.9
ACL TOP Family	Mean (% free PS)	CV% (Within run)	CV% (Total)
Normal Control	112.0	2.5	3.4
Special Test Control Level 1	64.9	1.9	3.2
Special Test Control Level 2	29.1	2.9	4.2
ACL TOP Family 50 Series ¹	Mean (% free PS)	CV% (Within run)	CV% (Total)
Normal Control	96.1	0.8	2.6
Special Test Control Level 1	59.8	1.0	2.4
Special Test Control Level 2	28.1	1.6	2.8

Correlation:				
System	slope	intercept	r	Comparative method
ACL Elite/Elite Pro	0.933	5.09	0.972	Free Protein S EIA
ACL TOP Family	1.000	-7.0	0.986	Free Protein S on ACL Advance
ACL TOP Family 50 Series ¹	0.940	0.63	0.997	Free Protein S on TOP 500 CTS

The precision and correlation results were obtained using specific lots of reagents and controls.

Detection limit:

System	
ACL Elite/Elite Pro	3.2 (% free Protein S)
ACL TOP Family/ACL TOP Family 50 Series ¹	7.3 (% free Protein S)

Linearity:

System	
ACL Elite/Elite Pro	10 - 135 (% free Protein S)
ACL TOP Family/ACL TOP Family 50 Series ¹	11 - 150 (% free Protein S)
If the linear range is exceeded samples can be manually diluted 1:2 with Factor diluent (200 µL of sample + 200 µL of Factor diluent) and reassayed. The printed results must be multiplied by 2 to return for the dilution.	
Instruments with Auto Rerun capability perform a 1:3 dilution and result correction automatically.	

† ACL TOP Family 50 Series = ACL TOP 350 CTS; ACL TOP 550 CTS; ACL TOP 750; ACL TOP 750 CTS; ACL TOP 750 LAS

Deutsch	Englisch	Franciais	Italiano
DEUTSCH - Packungsbeilage Version 04/2019	ENGLISH - Insert revision 04/2019	Français - Révision de la notice 04/2019	Italiano - Revisione del foglio 04/2019

Deutsch	Englisch	Franciais	Italiano
DEUTSCH - Packungsbeilage Version 04/2019	ENGLISH - Insert revision 04/2019	Français - Révision de la notice 04/2019	Italiano - Revisione del foglio 04/2019

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