Intended use

For the *in vitro* determination of Activated Partial Thromboplastin Time (APTT) in citrated plasma on IL Coagulation Systems as a general screening procedure for the evaluation of the intrinsic coagulation pathway and to monitor patients receiving heparin anticoagulant therapy.

Summary and principle

In the APTT test a contact activator is used to stimulate the production of Factor XIIa by providing a surface for the function of high molecular weight kininogen, kallikrein and Factor XIIa. This contact activation is allowed to proceed at 37°C for a specific period of time. Calcium is then added to trigger further reactions and the time required for clot formation is measured. Phospholipids are required to form complexes which activate Factor X and Prothrombin.

The APTT reagent included in the APTT-SP (liquid) kit contains synthetic phospholipids and silica to ensure a highly reproducible and stable product.^{1,2}

Prolonged clotting times may be observed in the following situations: deficiency of Factor XII, XI, X, IX, VIII, V, II, or fibrinogen, liver diseases, Vitamin K deficiency, presence of heparin, lupus anticoagulant or other inhibitor.^{3,4,5}

Composition

The APTT-SP kit consists of:

APTT reagent (Cat. No. 0020006310): 5 x 9 mLvials of colloidal silica dispersion with synthetic phospholipids, buffer and preservatives.

Calcium Chloride (Cat. No. 0019741900): 5 x 8 mL vials of calcium chloride (0.025 mol/L) with preservative.

PRECAUTIONS AND WARNINGS:

Hazard class: None

Hazard statements: None

Precautionary statements: None

Supplemental hazard information: None

Calcium Chloride contains sodium azide that may form explosive azides in metal plumbing. Use proper disposal procedures.

This product is For *in vitro* Diagnostic Use.

Preparation

APTT reagent: Shake silica dispersion vigorously for approximately 15 seconds or vortex for 5 seconds before use.

Calcium Chloride: The reagent is ready for use.

Reagent storage and stability

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.

APTT reagent: Opened reagent is stable 30 days at 2-8°C in the original vial or 5 days at 15°C on IL Coagulation Systems. No stirring is required. **Do not freeze.**

Calcium Chloride: Opened reagent is stable 30 days at 2-8°C in the original vial.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

Instrument/test procedures

Refer to the appropriate IL instrument's Operator's Manual and/or Application Manual for complete assay procedure instructions.

Specimen collection and preparation

Nine parts of freshly drawn venous blood are collected into one part 3.2% trisodium citrate. Refer to CLSI Document H21-A5 for further instructions on specimen collection, handling and storage.⁶

Additional reagents and control plasmas

The following are not supplied with the kit and must be purchased separately.

	Americas and Pacific Rim	Europe
	Cat. No.	Cat. No.
Calibration plasma	0020003700	0020003700
Normal Control	0020003120/0020003110	0020003110
Low Abnormal Control	0020003220/0020003210	0020003210
High Abnormal Control	0020003320/0020003310	0020003310
Factor diluent	0009757600	0009757600

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APTT-SP (liquid) - 0020006300

Quality control

Normal and abnormal controls are recommended for a complete quality control program.^{7,8} Normal Control, Low Abnormal Control and High Abnormal Control are designed for this program. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard *et al* for identification and resolution of out-of-control situations.⁹

Results

APTT results may be reported in seconds and/or ratios. Refer to the instrument's Operator's Manual for additional information.

Limitations/interfering substances

APTT results may be affected by many commonly administered drugs and further studies should be made to determine the source of unexpected abnormal results.

APTT results on IL Coagulation Systems are not affected by hemoglobin up to 100 mg/dL, triglycerides up to 1200 mg/dL and bilirubin up to 20 mg/dL.

Expected values

A normal range study was performed using APTT-SP reagent.

System	N	Range (units)
ACL 7000	63	24.3 - 35.0 (seconds)
ACL TOP® Family/ACL TOP Family 50 Series¹	120	25.4 - 36.9 (seconds)

Ranges were calculated as recommended by the International Federation of Clinical Chemistry (IFCC).¹⁰ These results were obtained using a specific lot of reagent.

Due to many variables which affect clotting times, each laboratory should verify its own normal range.

Performance characteristics

Precision:

Within run and total (run to run and day to day) precision was assessed over multiple runs using both normal and abnormal samples.

ACL 7000	Mean (APTT seconds)	CV % (Within run)	CV % (Total)
Normal Control	28.4	1.41	2.39
Low Abnormal Control	46.2	1.23	2.01
High Abnormal Control	72.2	1.88	2.80

ACL TOP Family	Mean (APTT seconds)	CV % (Within run)	CV % (Total)
Normal Control	29.0	1.1	2.5
Low Abnormal Control	43.6	2.4	3.6
High Abnormal Control	56.3	2.9	4.9

ACL TOP Family 50 Series¹	Mean (APTT seconds)	CV % (Within run)	CV % (Total)
Normal Control	29.6	0.6	0.9
Low Abnormal Control	41.4	0.8	1.4
High Abnormal Control	54.8	1.2	2.0

ACL Elite/Elite Pro	Mean (APTT seconds)	CV % (Within run)	CV % (Total)
Normal Control	29.3	1.1	2.0
Low Abnormal Control	46.1	1.9	4.1
High Abnormal Control	56.9	2.2	4.7

Correlation:

System	slope	intercept	r	Reference method*
ACL 7000	0.905	8.4	0.915	IL APTT-C
ACL 7000	0.827	6.1	0.942	IL APTT Lyophilized Silica
ACL TOP Family	0.934	4.3	0.987	APTT-SP on ACL Advance
ACL TOP Family 50 Series¹	0.956	1.6	0.999	APTT-SP on ACL TOP 500 CTS

These precision and correlation results were obtained using specific lots of reagents and controls.

*APTT Lyophilized Silica and APTT-C are not available in all countries.

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
Heparin Therapy¹¹


Therapeutic range: 0.2-0.4 U/mL

Heparin Dose-Response curve examples

The curve below is an example obtained using one lot of APTT-SP reagent and one Normal Plasma Pool spiked with graduated quantities of Bovine Unfractionated Heparin.

The curve below is an example obtained using one lot of APTT-SP reagent and plasma samples from patients undergoing heparin therapy (Sodium Unfractionated heparin), to show the correspondence, for each level of heparin, between the heparin anti Xa activity and the relevant APTT response.

Due to many variables (i.e. different source of heparin) which may affect the clotting times, each laboratory should establish its own heparin therapeutic range.¹²

For further indications on therapeutic ranges and duration of treatment make reference to local guidelines.

Intrinsic Factor Sensitivity

Studies have shown APTT-SP to be sensitive to decreased concentration of intrinsic factors resulting in an abnormal APTT value when factors VIII, IX, XI and XII levels were in the 30-60% range.

† ACL TOP Family 50 Series = ACL TOP 350 CTS; ACL TOP 550 CTS; ACL TOP 750; ACL TOP 750 CTS; ACL TOP 750 LAS

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