Tina-quant Lipoprotein(a) Gen.2

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Letteratura

- 1 Siekmeier R, Schamagl H, Kostner GM, et al. Lipoprotein(a) - Structure, Epidemiology, Function and Diagnostics of a Cardiovascular Risk Marker. *The Open Clin Chem J* 2008;1:79-91.
- 2 Kamstrup PR. Lipoprotein(a) and Ischemic Heart Disease- A Causal Association? A review. *Atherosclerosis* 2010 Jul;211(1):15-23.
- 3 Genser B, Dias KC, Siekmeier R, et al. Lipoprotein(a) and Risk of Cardiovascular Disease - A Systematic Review and Meta Analysis of Prospective Studies. *Clin Lab* 2011;57(3-4):143-156.
- 4 Kamstrup PR, Tybjaerg-Hansen A, Steffensen R, et al. Genetically Elevated Lipoprotein(a) and Increased Risk of Myocardial Infarction. *JAMA* 2009;301(22):2331-2339.
- 5 Clarke R, Peden JF, Hopewell JC, et al. Genetic Variants Associated with Lp(a) Lipoprotein Level and Coronary Disease. *N Engl J Med* 2009 Dec;361(26):2518-2528.
- 6 Dati F, Tate JR, Marcovina SM, et al. First WHO/IFCC International reference Reagent for Lipoprotein(a) for immunoassay - Lp(a) SRM2B. *Clin Chem Lab Med* 2004;42(6):670-676.
- 7 Marcovina SM, Albers JJ, Gabel B, et al. Effect of the Number of Apolipoprotein (a) Kringle 4 Domains on Immunochemical Measurements of Lipoprotein(a). *Clin Chem* 1995 Feb;41(2):246-255.
- 8 Marcovina SM, Albers JJ, Wijsman E, et al. Differences in Lp(a) Concentrations and Apo(a) Polymorphs Between Black and White Americans. *J Lipid Res* 1996 Dec;37(12):2569-2585.
- 9 Reiner Ž, Catapano AL, De Backer G, et al. ESC/EAS Guidelines for the management of dyslipidaemias. *Eur Heart J* 2011;32:1769-1818.
- 10 Simó JM, Camps J, Gómez F, et al. Evaluation of a Fully Automated Particle-enhanced Turbidimetric Immunoassay for the Measurement of Plasma Lipoprotein(a). Population-Based Reference Values in an Area with Low Incidence of Cardiovascular Disease. *Clin Biochem* 2003 Mar;36(2):129-134.
- 11 National Committee for Clinical Laboratory Standards, Procedures for the handling and Processing of Blood Specimens, Approved Guideline, NCCLS publication H18-A, Villanova, 1990.
- 12 Simó JM, Camps J, Vilella E, et al. Instability of Lipoprotein (a) in Plasma Stored at -70 °C: Effects of Concentration, Apolipoprotein (a) Genotype, and Donor Cardiovascular Disease. *Clin Chem* 2001 Sep;47(9):1673-1678.
- 13 Sgoutas DS, Tuten T. Effect of Freezing and Thawing of Serum on the Immunoassay of Lipoprotein(a). *Clin Chem* 1992;38(9):1873-1877.
- 14 Marcovina SM, Albers JJ, Scanu AM, et al. Use of a reference Material Proposed by the International federation of Clinical Chemistry and Laboratory Medicine to Evaluate Analytical methods for the Determination of Plasma Lipoprotein (a). *Clin Chem* 2000 Dec;46(12):1956-1967.
- 15 Anne Langsted, Pia R. Kamstrup, Borge G. Nordestgaard. High lipoprotein(a) and high risk of mortality. *European Heart Journal* (2019) 40:2760-2770.
- 16 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- 17 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". *Eur J Clin Chem Clin Biochem* 1996;34:385-386.
- 18 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.
- 19 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *Clin Chem Lab Med* 2007;45(9):1240-1243.
- 20 Marcovina SM, Koschinsky ML. A Critical Evaluation of the Role of Lp(a) in Cardiovascular Disease: Can Lp(a) Be Useful in Risk Assessment? *Semin Vasc Med* 2002 Aug;2(3):335-344.
- 21 Shai I, Rimm EB, Hankinson SE, et al. Lipoprotein (a) and Coronary Heart Disease Among Women: Beyond a Cholesterol Carrier? *Eur Heart J* 2005;26:1633-1639.
- 22 Nordestgaard BG, Chapman MJ, Ray K, et al. Lipoprotein (a) as a cardiovascular risk factor: current status. *Eur Heart J* 2010 Dec;31(23):2844-2853.
- 23 Marcovina SM, Koschinsky ML, Albers JJ, et al. Report of the National Heart, Lung, and Blood Institute Workshop on Lipoprotein (a) and Cardiovascular Disease: Recent Advances and Future Directions. *Clin Chem* 2003 Nov;49(11):1785-1796.
- 24 Tsimikas S, Clopton P, Brilakis ES, et al. Relationship of Oxidized Phospholipids on Apolipoprotein B-100 Particles to Race/Ethnicity, Apolipoprotein (a) Isoform Size, and Cardiovascular Risk Factors: Results From the Dallas Heart Study, *Circulation* 2009 Apr;119(13):1711-1719.
- 25 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

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© 2019, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim