LabUSticks 11F 12F 14F

ENGLISH

LabUSticks 11F/12F/14F Test Strips for Urine

READ THIS PACKAGE INSERT CAREFULLY BEFORE USE

For In Vitro Diagnostic Use Only

INTENDED USE

LabUSticks 11F/12F/14F Test Strips for Urine provide tests for the semi-quantitative measurement of ascorbic acid, nitrite, microalbumin, leukocytes, creatinine, ketones, urobilinogen, bilirubin, glucose, proteins, specific gravity, pH, blood and calcium in urine (refer to the carton or bottle label to see which tests are included on the strip you are using). Use only with the UNAMAX automatic urine analyzer and the UNAMAX Lite urine analyzer. LabUSticks 11F/12F/14F Test Strips for Urine are intended for in vitro diagnostic use by healthcare professionals.

SUMMARY

LabUSticks 11F/12F/14F Test Strips for Urine consist of a plastic strip equipped with reagent paper pads and a calibration pad. This feature facilitates the measurement of multiple urine constituents and use for everyday diagnosis and group examinations. The calibration pad, which is not impregnated with reagents, allows for automatic instrumental correction of deviations from the natural urine color in order to obtain accurate results.

TEST PRINCIPLES AND LIMITATIONS

Ascorbic acid: The test involves the decolorization of Tillman's reagent. A false positive reaction may occur in the presence of other reducing agents.

Nitrite: This test is based on the principle of the Griess test and is specific for nitrites. The development of any degree of uniform pink color should be interpreted as a positive. The nitrite test suggests the presence of 10⁶ or more organisms per mL, but color development is not proportional to the number of bacteria present. A negative result does not in itself prove that there is no significant bacteria. Negative results may occur when urinary tract infections are caused by organisms which do not contain reductase to convert a nitrate to a nitrite, even when urine has been retained long enough (4-8 hours) for nitrate reduction to occur, or when dietary nitrate is absent, even if organisms containing reductase are present and bladder incubation is ample. Ascorbic acid concentrations of 1.4 mmol/L or greater may cause false negative results with specimens containing nitrite ion concentrations of 43 µmol/L or less. Owing to the existence of nitrates, a positive result may not indicate a urinary system infection.

Microalbumin: The reaction of albumin is more sensitive than the reaction of globulin, hemoglobin, Bence-Jones protein, and mucin, thus a negative result does not rule out the existence of the above-mentioned proteins in urine. When the results are 20-200 µg/L, microalbuminuria is indicated, and when the results are beyond 200 µg/L, clinical albuminuria is indicated. This action is usually affected by creatinine and hemoglobin, etc.. Alkaline urine may cause false positive results.

Leukocytes: This test reveals the presence of granulocyte esterase. These esterase cleave an indoxyl ester, and the indoxyl liberated reacts with a diazonium salt to produce a violet dye. Leukocyte esterase results may be positive in the absence of observable cells if the leukocytes have lysed. Positive results may occasionally be found with random specimens from females due to contamination of the specimen by vaginal discharge. Elevated glucose concentrations (>55 mmol/L) or high specific gravity may reduce the test results. The presence of cephalopins, cephalothin, tetracycline may cause decreased reactivity, and high levels of the drug may cause a false negative reaction. The test area does not react with lymphocytes. Reactivity may also vary with temperature.

Creatinine: This test is based on the principle of a displacement reaction. Creatinine displaces the dyes from the compound of metallic chloride and acid dyes, causing a color change from green to yellow. Daily creatinine excretion, related to muscle mass of the human body, is usually constant. Some compounds, physical properties and high concentration yellow pigments may affect the test results.

Ketones: This test is based on the principle of Legal's test and is more sensitive to acetoacetic acid than acetone. The reagent area does not react with β-hydroxybutyric acid. Some high specific gravity or low pH urine may give reactions up to and including trace levels. Normal urine specimens usually yield negative results with this reagent. False positive results (Trace) may occur with highly concentrated urine specimens containing large amounts of urea or levodopa. Urine specimens containing a yellow or reddish discoloration of the test paper that can interfere with interpreting the result. Ascorbic acid concentrations of 1.4 mmol/L or greater may cause false negatives. **Glucose:** This test, which is based on the specific glucose oxidase/peroxidase reaction, is specific for glucose. No substance excreted in urine other than glucose is known to give a positive result. Ascorbic acid of more than 1.4 mmol/L and/or high ketone concentrations (8 mmol/L) may cause false negatives for specimens containing small amounts of glucose (5.5 mmol/L). The reactivity of the glucose test decreases as the specific gravity of the urine increases. False positive reactions may be caused by hypochlorite or peroxide (bleaching agents). Reactivity may also vary with temperature.

Proteins: This test is based on the principle of the protein error of pH indicator. The reagent area is more sensitive to albumin. An elevated pH (pH 9) will not affect the test. The residues of disinfectants containing quaternary ammonium groups or chlorhexidine present in the urine vessel may lead to false positive results. **Specific gravity:** This test contains a detergent and Bromothymol blue that indicates the presence of cationic constituents in the urine by changing color from green to yellow. The specific gravity test allows for determining the urine specific gravity between 1.005 and 1.030. In general, it correlates within 0.005 with values obtained with the refractive index method. Strips are automatically adjusted for pH by the instrument when pH 7.0 or pH 5.0. Highly buffered alkaline urine may cause low readings compared to other methods. Elevated specific gravity readings may be obtained in the presence of moderate quantities of protein (5 g/L). **pH:** This test contains a mixed indicator which generates an evident color change between pH 5.0 and pH 9.0.

Blood: Hemoglobin and myoglobin catalyze the oxidation of the indicator by means of organic hydroquinone contained in the test paper. This test is highly sensitive to hemoglobin and thus complements the microscopic examination. The sensitivity of this test may be reduced in urine with high specific gravity. The test is equally sensitive to both myoglobin and hemoglobin (hemoglobin concentration of 150-620 µg/L is approximately equivalent to 5-15 intact red blood cells per microliter). Captopril and Etodolac (Lodine) may also cause decreased reactivity. Blood is often found in the urine of menstruating females. Certain oxidizing contaminants, such as hypochlorite, may produce false positive results. Microbial peroxidase associated with urinary tract infections may cause a false positive reaction. Ascorbic acid concentrations of 1.4 mmol/L or greater may cause false negatives at the trace levels.

Calcium: This test is based on the calcium ion in urine which reacts with oCPC to produce a color change. A high concentration of magnesium ion in urine may affect the test. **Microalbumin-to-creatinine (A:C) ratio:** The ratio obtained by combining microalbumin and creatinine can reduce the stochastic error. Normally, the A:C ratio is lower than 30 mg/g (3.4 mg/mmol). When the ratio is between 30 mg/g (3.4 mg/mmol) and 300 mg/g (33.9 mg/mmol), it is abnormal, and when it is higher than 300 mg/g (33.9 mg/mmol), it is severely abnormal.

TEST PROCEDURE

- Additional materials required: UNAMAX automatic urine analyzer or UNAMAX Lite urine analyzer. Detailed operating instructions can be found in the instruction manual of the instrument.
- When using UNAMAX Lite urine analyzer, please ensure that the urine test strip is fully dipped into the urine sample or use a dropper to drip the urine sample evenly on the test strip. When using UNAMAX automatic urine analyzer, please ensure that the urine sample is dropped against the rim of the urine container to remove excess urine. Blot the longer edge of the strip with absorbent paper. Avoid any runner contamination from adjacent reagent pads.
- When interpreting the result with the urine analyzer, please follow the instructions in the relative operating manual.

PRECAUTIONS

- The test is intended for use by health care professionals.
- The color chart on the label is intended as a reference only.
- Collect fresh urine specimens in clean, dry containers. Do not expose urine specimens to sunlight as this induces oxidation of the bilirubin and urobilinogen and hence gives artificially lower results for these two parameters.
- Do not touch the reagent pads. The strips should be kept clean to avoid contamination.
- Do not remove desiccants. Immediately rewrap tightly after removing the reagent strip: lengthy exposure (more than 5 minutes) to air humidity can easily give rise to inaccurate test results.
- Do not use reagent strips after their expiry date, and do not use deteriorated, discolored or blackened test strips.
- If refrigerated, return strips to room temperature before use.
- False positive reactions for blood and glucose may result from residues of strongly oxidizing disinfectants in the specimen collection vessel. Do not add preservatives to the urine. Avoid contamination by volatile chemicals.
- Used test strips must not be reused and should be disposed of as general medical waste.
- As a rule, diagnosis should not be based on the results of a single test established in the context of all other medical findings. Full knowledge of the effects of drugs or their metabolites on the individual tests is not yet available. In doubtful cases, it is therefore advisable to repeat the test after discontinuing a specific drug. Large amounts of ascorbic acid in the urine can produce artificially low to false-negative results for glucose, blood, nitrites, and bilirubin.

STORAGE AND STABILITY

Store at temperatures between 2 °C and 30 °C. Store only in the original bottle, avoiding humidity, direct sunlight, or heat. Unused strips that remain in the original capped container are stable for 5 months after first opening. Do not use reagent strips after the expiry date printed on the vial label. If unopened the product is valid for 18 months.

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